



**LUNG
SASK**

RESPICLICK™
Inhaler Delivery Device

HOW TO USE A RESPICLICK™

1. Hold the inhaler upright.
2. Pull the cap downwards until you hear a click.
3. Breathe out, away from the mouthpiece.
4. Seal your lips around the mouthpiece.
5. Tilt your chin up slightly and quickly breathe in until your lungs are full.
6. Remove the inhaler from your mouth and hold your breath for 5-10 seconds
7. Breathe out, away from the mouthpiece.
8. Close the cap.
9. Rinse your mouth with water and spit the water out.




1. Hold the inhaler upright



2. Pull cap downwards until you hear a click



3. Breathe out



4. Seal lips around mouthpiece




5. Tilt chin up slightly and breathe in quickly



6. Hold breath (5-10 secs)



7. Breathe out
8. Close the cap



9. Rinse your mouth with water

CARE FOR A RESPICLICK™

1. Discard the inhaler 30 days after removing from the foil pouch or when the counter reads “0”.
2. Once the cap is open, do not close until the dose of medicine has been inhaled, or this will waste a dose.
3. Keep the inhaler dry. If you need to clean it, gently wipe the mouthpiece with a damp cloth or dry tissue.
4. When the dose counter reaches 20, the number of the colors will change to red. This is a reminder to refill your prescription.
5. The RespiClick™ counts down by twos (i.e., 60, 58, 56, etc.). When the counter shows “0”, the background will change to red and the “0” will be a black font. This means the inhaler is empty.
6. The RespiClick™ delivers your dose of medicine as a very fine powder that you may or may not taste or feel. Do not take an extra dose from the device even if you do not taste or feel the medicine.
7. Always check the inhaler’s expiration date.

Important information:

RespiClick™ is the name of the inhaler delivery device.
Medication in this device has a brand name and generic name.

Inhalers available in the RespiClick™:

Aermony

Generic name: fluticasone propionate

Dose per puff: 55, 113, or 232 mcg

Duration: 12 hours



For video instruction:
lungsask.ca/inhaler-resources



lungsask.ca

